



BRIEF INTRODUCTION

Zhejiang Bangli Medical Products Co., Ltd.(registered fund \$1.3 million, annual sales amount \$30 million and keep growing 40% every year)established in 2003 ,is a foreign market-oriented manufacturer of first-aid bandages, surgical plaster, cooling gel sheets, nasal strip, Pain relief patch , cold & hot patch, wound Dressing plaster and other related medical device products.

At present, we can produce 4,000,000pcs of first aid bandages, 200,000rolls of cohesive bandage, 200,000rolls of surgical tape, 200,000pcs of cooling gel patches, 200,000pcs of nasal strips, 300,000pcs of wound dressing per day and our production capacity is enlarging step by step. We have exported to European and American countries and regions such as the USA,UK,Germany,Swiss,Italy and we also have enlarged to the Middle East, SOUTH EAST ASIA, SOUTH AFRICA AND other countries and regions. We welcome all the customers in every corner of the world.

To ensure quality, we have built our own lab, sterilization workshop and one hundred thousand clean room workshop. We have CE and FDA certificates for our products. In addition, we have developed a practical quality management system and approved by ISO 9001:2000 and ISO 13485 .

Rooting on the culture cooperation of“Harmony, Share and Innovation”Zhejiang Bangli Medical Products Co.,Ltd has paid much attention to the talents and product innovation, Aspiring of the spirit of serving the society, we will do every endeavor to the quality and service for our every customer at all the time .

BL 邦立
H/E/A/L/T/H / C/A/R/E



COMPANY PROFILE



Section of production line



Section of packing room



Arab Health in Dubai



Foreign trade department



Sample display room




Medica In Germany

has passed ISO9001 Quality System Certification.
ISO13485 Quality Certification of Germany TUV Authentication Institution Organization.
GMP Quality System Certification. CE Certification and the FDA Certification.
It is owned self-support import and export rights.



企业组织结构

ENTERPRISE ORGANIZATION STRUCTERE



CATALOGUE



01

- 03···FIRST AID BANDAGE
- 06···FEVER RELIEF PATCH
- 07···COLD & HOT THERAPY
- 08···PAIN RELIEF PATCH

03

- 23···COHESIVE BANDAGE
- 24···ADHESIVE BANDAGE
- 25···SPORTS TAPE
- 28···MOSQUITO REPELLENT PATCH
- 28···NASAL INHALANT
- 29···CORN PLASTER
- 30···ESSENTIAL BALM
- 30···ESSENTIAL OIL
- 31···FOAM DRESSING
- 32···SILICONE PLASTER ON ROLL
- 33···ONE-PIECE COLOSTOMY BAG
- 34···TWO-PIECE COLOSTOMY BAG
- 34···THE LAYERS CAN BE CUSTOMIZED

02

- 12···SURGICAL TAPE
- 16···NASAL STRIPS
- 16···VAPOUR PATCH
- 17···WOUND DRESSING
- 19···EXTENSIBLE PLASTER ROLL
- 19···EYE PATCH
- 20···FUNCTIONAL WOUND DRESSING
- 22···ACNE PLASTERS
- 22···COLOSTOMY CARE

ZHEJIANG BANGGLI
MEDICAL
PRODUCTS CO.,LTD
PROFESSIONAL FACTORY
WITH 20 YEARS' EXPERIENCE



FIRST AID

01

- 03···FIRST AID BANDAGE
- 06···FEVER RELIEF PATCH
- 07···COLD & HOT THERAPY
- 08···PAIN RELIEF PATCH

FIRST AID



FIRST AID BANDAGE

FLAT FABRIC BANDAGE



Item No.	Specification	Pcs/box
BL-001	72mmx19mm	10pc/small box 30small boxes/display box
BL-002	72mmx19mm	100
BL-003	70mmx18mm	100

RED COTTON ADHESIVE STRIP



Item No.	Specification	Pcs/box	Boxes/ctn
BL-004	60mmx15mm	100	100
BL-005	76mmx38mm	20pcs/size	100
	76mmx25mm		
	72mmx19mm		
	76mmx45mm		
BL-006	38mmx38mm	100	100
	65mmx19mm		

BLUE ELASTIC FABRIC BAND-AID , WITH DETECTIVE METAL PAD



USED FOR FOOD AND DRUG INDUSTRY

Item No.	Specification	Pcs/box	Boxes/ctn
BL-007	50mmx45mm	40	50
BL-007-L	76mmx45mm	25	50
BL-008	76mmx38mm	40	50
BL-009	76mmx25mm	100	50
BL-010	76mmx22mm	50	50
BL-011	76mmx50mm	25	50

LONG STRIPS



PE , NON-WOVEN , ELASTIC FABRIC MATERIAL

Item No.	Specification	Pcs/box	Boxes/ctn
BL-012	6cm x 1m	10	50
BL-013	7cm x 1m	10	50
BL-014	8cm x 1m	10	50

HIGH ELASTIC FABRIC WOUND PLASTER




Item No.	Specification	Pcs/box	Boxes/ctn
BL-015	76mmx22mm	50	50
BL-016	76mmx45mm	25	50
BL-017	76mmx38mm	40	50
BL-018	76mmx50mm	25	50
BL-019	50mmx45mm	40	50

HEAVY WEIGHT FABRIC BAND-AID



Item No.	Specification	Pcs/box	Boxes/ctn
BL-020	76mmx22mm	50	50
BL-021	76mmx45mm	25	50
BL-021-S	50mmx45mm	40	50
BL-022	76mmx38mm	40	50
BL-023	38mmx22mm	50	100
BL-024	76mmx50mm	25	50

NON-WOVEN ADHESIVE PLASTER



WHITE COLOR / WITH HYPOALLERGENIC GLUE

Item No.	Specification	Pcs/box	Boxes/ctn
BL-025	76mmx19mm	100	100
BL-026	76mmx19mm	30	480

PE, PVC ,MATERIAL BANDAGE



WE HAVE SKIN AND TRANSPARENT COLOR

Item No.	Specification	Pcs/box
BL-027	Skin color 76mmx19mm	100
BL-028	Skin color Assorted pack	100
BL-029	Skin Color Assorted 30	30
BL-030	Transparent 72mmx19mm	100
BL-031	Transparent Assorted pack	100
BL-032	Transparent 72mmx19mm	30
BL-033	Transparent Assorted pack	30



PU BANDAGE



Item No.	size	Pcs/box
BL-034	Transparent PU	100pcs
	38mmx38mm	
	assorted size 65mmx25mm	
	76x19mm	

FUNKIDS CARTOON BANDAGE , PRINTED PE OR PVC .



Item No.	size	Pcs/box
BL-035-TC	72mmx19mmwith printed, Tin box pack	25pcs
BL-036-C	76mmx19mmwith printed, Papercard box pack	25pcs
BL-037-C	76mmx19mm with printed.Papercard box pack	40pcs
BL-038-PC	72mmx19mmwith printed. Plastic box pack	5pcs

BL-039 BUTTERFLY CLOSURES

Medium 45mmx10mm 100pcs/box
Large 68mmx12mm 100pcs/box




BL-040 FACE STICKER

Good choice for promotion and sports

Specialize in printing high quality and fast product whether you need band stickers, school stickers ,or stickers for your church, club,business or organization.

30mmx45mm
38mmx38mm



FEVER RELIEF PATCH

COOLING GEL PATCH

Material: Non-woven fabric with hydro-gel and opp film

Areas of use: Effective to the fever and relieve pain.

Size: 4cmx10cm 5cmx12cm 5cmx11cm 10cmx14cm

Smell: Apple. Orange. Menthol. Candy Smell etc.

Color: All kinds of colors as you wish



BL-041 PATCH FOR KIDS



BL-042 PATCH FOR ADULT



BL-043 PATCH FOR ADULT



BL-044 PATCH FOR BACK

Item No.	size	patches/box	patches/carton	carton size (cm)	G.W(kgs)
BL-041	4cmx10cm	6	300	52x38x13	5
BL-042	5cmx11cm	6	300	52x38x13	6
BL-043	5cmx12cm	6	300	52x38x13	8
BL-044	10cmx14cm	6	180	46.5x41.5x7	9



COLD & HOT THERAPY

COLD&HOT PATCH

Material: Elastic non-woven and Hydrophilic polymer gel and release layer

Areas of use:

Pain on arthritis, backache, tendonitis, muscle strains, sprains, bruises and cramps




Item No.	size	patches/box	patches/carton	carton size (cm)	G.W(kgs)
BL-045	8cmx12cm	5	120	27x37x13.5	4
BL-046	10cmx14cm	5	120	32x35x13.5	5
BL-047	10cmx20cm	5	120	32x35x13.5	7

PAIN RELIEF PATCH

CAPSICUM PLASTER


Composition: Fabric + capsicum, menthol, crystals, Ginger and borneol etc.

Use: Relieve rheumatism & muscle pain like lumbago, chest, headache, toothache.



TO ALLOW FREE PASSAGE OF AIR.
TO IMPROVE ADHESION TO THE SKIN.
WHITE, SKIN, GRAY, ORANGE OR ANY COLOR

Item No.	size	pcs/box	pcs/ctn
BL-049	4cmx6.4cm	200	8000
BL-050	10cmx18cm	20	800
BL-051	12cmx18cm	24	960





MEDICATED PLASTER

Composition: Fabric(cotton), mentol, crystals Ginger and borneol etc.

Use: Dispel rheumatism, muscular pain like lumbago, chest and lung.




MORE EFFECTIVE AFTER A HOT BATH.

DO NOT APPLY TO IRRITATED SKIN.

WHITE COLOR, SKIN, GRAY, ORANGE, OR ANY COLOR U NEED.

Item No.	size	pcs/box	pcs/ctn
BL-052	4cmx6.4cm	200	8000
BL-053	10cmx18cm	20	800
BL-054	12cmx18cm	24	960



ZHEJIANG BANGLI
MEDICAL
PRODUCTS CO., LTD
PROFESSIONAL FACTORY
WITH 20 YEARS' EXPERIENCE

MEDICAL
TAPES & DRESSING

MEDICAL TAPES & DRESSING

02

- 12···SURGICAL TAPE
- 16···NASAL STRIPS
- 16···VAPOUR PATCH
- 17···WOUND DRESSING
- 19···EXTENSIBLE PLASTER ROLL
- 19···EYE PATCH
- 20···FUNCTIONAL WOUND DRESSING
- 22···ACNE PLASTERS
- 22···COLOSTOMY CARE



SURGICAL TAPE

ZINC OXIDE PLASTER

Description

- * Material: 100% cotton with hypoallergenic glue
- * Good stickness, low allergic, easy tear
- * High tensile strength, provide strong and controlled support
- * Good weatherability * Size as customized



BL-057 **
SIMPLE PACKING



BL-058 **
PLASTIC WING



BL-059-1



BL-059-2



BL-060 **
METAL TIN PACKING



BL-S **
PERFORATED

** SPOOL PACKAGED, DUST-PROOF

Pack	simple pack	Plastic spool pack	Metal tin pack
1.25cmx5m	24rolls/box, 30boxes/ctn	24rolls/box, 12boxes/ctn	24rolls/box, 12boxes/ctn
2.5cmx5m	12rolls/box, 30boxes/ctn	12rolls/box, 30boxes/ctn	12rolls/box, 30boxes/ctn
5cmx5m	6rolls/box, 30boxes/ctn	6rolls/box, 30boxes/ctn	6rolls/box, 30boxes/ctn
7.5cmx5m	6rolls/box, 30boxes/ctn	6rolls/box, 30boxes/ctn	6rolls/box, 30boxes/ctn
10cmx5m	6rolls/box, 30boxes/ctn	6rolls/box, 30boxes/ctn	6rolls/box, 30boxes/ctn
CAN BE ANY LENGTH , SUCH AS: 5YARDS, 10 YARDS, 10 M AND SO ON			

SILK TAPE

Description

- * Material: Silk cloth with hypoallergenic glue
- * Very high tack, Good stickness
- * Silk Plaster with/without zigzag



BL-061-A *** PACK IN TRANSPARENT
PLASTIC COVER AND PUT IN INDIVIDUAL BOX



BL-061 *** SIMPLE PACKING
(CAN BE PLAIN EDGE OR ZIGZAG)



BL-062 *** PLASTIC SPOOL PACKAGED, DUST-PROOF

Pack	simple pack	Plastic spool pack
1.25cmx5m	24rolls/box, 30boxes/ctn	1 roll/Plastic spool, 1 roll/box, 480boxes/ctn
2.5cmx5m	12rolls/box, 30boxes/ctn	1 roll/Plastic spool, 1 roll/box, 360boxes/ctn
2.5cmx9.1m	12rolls/box, 30boxes/ctn	1 roll/Plastic spool, 1 roll/box, 360boxes/ctn
5cmx5m	6rolls/box, 30boxes/ctn	1 roll/Plastic spool, 1 roll/box, 180boxes/ctn
5cmx9.1m	6rolls/box, 30boxes/ctn	1 roll/Plastic spool, 1 roll/box, 180boxes/ctn
7.5cmx5m	6rolls/box, 30boxes/ctn	1 roll/Plastic spool, 1 roll/box, 180boxes/ctn
10cmx5m	6rolls/box, 30boxes/ctn	1 roll/Plastic spool, 1 roll/box, 120boxes/ctn
CAN BE ANY LENGTH , SUCH AS: 5YARDS, 10 YARDS, 10 M AND SO ON		



PE TAPE

Description

- * Waterproof
- * Breathable
- * Good skin adhesion
- * Transparent and blue color



Transparent



Semi-transparent

BL-063 *** SIMPLE PACKING



BL-064 ***

PLASTIC SPOOL PACKAGED,DUST-PROOF



BL-065 ***

PE WITH CUTTER DISPENSER

Pack	simple pack	Plastic spool pack	with cutter dispenser
1.25cmx5m	24rolls/box,30boxes/ctn	24rolls/box,12boxes/ctn	1roll/cutter,24cutter/box,30boxes/ctn
2.5 cmx5m	12rolls/box,30boxes/ctn	12rolls/box,30boxes/ctn	1roll/cutter,12cutter/box,30boxes/ctn
5cmx5m	6rolls/box,30boxes/ctn	6rolls/box,30boxes/ctn	1roll/cutter,6cutter/box,30boxes/ctn
7.5cmx5m	6rolls/box,30boxes/ctn	6rolls/box,30boxes/ctn	1roll/cutter,6cutter/box,30boxes/ctn
10cmx5m	6rolls/box,30boxes/ctn	6rolls/box,30boxes/ctn	1roll/cutter,6cutter/box,30boxes/ctn
CAN BE ANY LENGTH , SUCH AS: 5YARDS, 10 YARDS, 10 M ETC.			

NON-WOVEN PAPER TAPE

Description

- * Material:Non-woven with hypoallergenic glue
- * Breathable soft and comfortable
- * High tensile strength,non-slip
- * Remove easily,non-sticky to skin or hair,no degumming,leaves no residue



BL-066 *** SIMPLE PACKING



BL-067 *** WITH CUTTER DISPENSER PACKING

Pack	simple pack	with cutter dispenser
1.25cmx5m	24rolls/box,30boxes/ctn	1roll/cutter,24cutter/box,30boxes/ctn
2.5 cmx5m	12rolls/box,30boxes/ctn	1roll/cutter,12cutter/box,30boxes/ctn
5cmx5m	6rolls/box,30boxes/ctn	1roll/cutter,6cutter/box,30boxes/ctn
7.5cmx5m	6rolls/box,30boxes/ctn	1roll/cutter,6cutter/box,30boxes/ctn
10cmx5m	6rolls/box,30boxes/ctn	1roll/cutter,6cutter/box,30boxes/ctn
CAN BE ANY LENGTH , SUCH AS: 5YARDS, 10 YARDS, 10 M ETC.		



NASAL STRIPS

Expand nasal passages and improves airflow by 31%;

Normally no rash or skin irritation

Material: PE/Non-woven

Pack: 20pcs or 50pcs/box. 240 boxes/ctn

Size: 55mmx18mm 58mmx18mm 66mmx19mm 73mmx19mm



BL-068



BL-069-1
PE MATERIAL




BL-069-2



BL-070



BL-071-1



VAPOUR PATCH

- Cool easy Breathing
- Slow release patch over 8 hours
- Clear and Easy Breathing



BL-072



WOUND DRESSING

COMMON DRESSING

package 1pc/envelope 50pcs/box, 20boxes/ctn

Size: 5cmx7cm 8cmx10cm 9cmx10cm 10cmx10cm

15cmx10cm 20cmx10cm 25cmx10cm 30cmx10cm 35cmx10cm etc

NON-WOVEN WOUND DRESSING PAD



BL-073



PU WOUND DRESSING PATCH



BL-074



PUNCTURE FIXING IV DRESSING

Size: 6cmx8cm

Application: For fixing Y or Butterfly catheter or catheter with drug kettle


Pack: 50pcs/box 40 boxes/ctn



Spunlace nonwoven fabric coated with acrylic adhesive can stick on the skin firmly, it will not come off from the skin when meeting with water, it is hypoallergenic and easy to be peeled off.

"U" cutting style enhances the fixing of the catheter.

The pad has good absorption, can absorb the exudation from the pinhole well.



Item No.	Size	Packing
BL-075	4.4cmX6 cm	100pcs×20Box 2000pcs/ctn
	6cm×7cm	100pcs×20Box 2000pcs/ctn
	10m×12cm	50pcs/box,24boxes/ctn
	5cm×5.7cm	100pcs×20Box 2000pcs/ctn
	6.5cmX7cm	100pcs/box, 24boxes/ctn
	6cm×7cm	100pcs×20Box 2000pcs/ctn
	7cm×9cm	100pcs×20Box 2000pcs/ctn
	6cm×8cm	100pcs/box, 24boxes/ctn
	7cm×9.5cm	100pcs/box, 24boxes/ctn
	8.5cmX10.5cm	100pcs/box, 24boxes/ctn
	7cm×8.5cm	100pcs/box, 24boxes/ctn

EXTENSIBLE PLASTER ROLL

Non-woven

PU



BL-076

BL-076-1

Size: 2.5cmx10m, 5cmx10m, 10cmx10m, or customized.

EYE PATCH BL-077

Non-adhesive island dressing with shaped pad to allow eyelid movement;

Easy to apply. Lightweight, hypoallergenic and breathable; Gentle adhesive edge;

Comfortably behind glasses;



With design


Size: 70mmx54mm; 80mmx57mm;95mmx65mm

Color: White, Red, Yellow or customized



FUNCTIONAL WOUND DRESSING

Hydrocolloid dressing is a kind of elastic dressing made of medical hot-melt adhesive with CMC hydrophilic granules.



Item No.	Size	Packing
BL-079	44mmx69mm	1pcs/bag, 10bags/box,40boxes/ctn
	55mmx37mm	1pcs/bag, 10bags/box,40boxes/ctn
	20mmx60mm	1pcs/bag, 10bags/box,40boxes/ctn
	20mmx65mm	1pcs/bag, 10bags/box,40boxes/ctn
	28mmx69mm	1pcs/bag, 10bags/box,40boxes/ctn
	42mmx62mm	1pcs/bag, 10bags/box,40boxes/ctn
	28mmx46mm	1pcs/bag, 10bags/box,40boxes/ctn



Characteritics:

1. excellent absorption capability
2. keep the wound a moist environment, accelerate wound healing process, effectively reduce pains and wound care frequency
3. water-proof, air permit, keep the wound from outer contamination
4. self-adhesive and with suitable elastic, comfortable to wear
5. please change the dresssing when it turns white
6. easy to move off, and no adhesive with wound, so effectiely reduce the possible damage and pains.

HYDROCOLLOID DRESSING PAD



Item No.	size	Packing
BL-078	5cmx5cm	50pcs/box, 1800pcs/ctn
	7.5cmx7.5cm	25pcs/box, 900pcs/ctn
	10cmx10cm	25pcs/box, 500pcs/ctn
	15cmx15cm	10pcs/box, 480pcs/ctn
	20cmx20cm	10pcs/box, 300pcs/ctn
	8cmx12cm	25pcs/box, 900pcs/ctn
	15cmx18cm	25pcs/box, 300pcs/ctn

HYDROCOLLOID PLASTER/BLISTER PLASTER

Manage wound exudate;
Create a moist wound environment;
Deliver long lasting antimicrobial protection for up to seven days.

Size: 44mmx69mm, 55mmx37mm, 20mmx60mm, 20mmx65mm, 69mmx28mm, 55mmx37mm, 44mmx69mm, 50mmx50mm, 30mmx60mm, 25mmx65mm, 15mmx51mm, 46mmx28mm etc.



BL-079

FRON EASE


BL-080



Size: 26mmx15mm



ACNE PLASTERS



Item No.
BL94



Description
24grain



Packing
1pcs/pouch,
50 pouches/bag;
40bags/ctn

ZHEJIANG BANGLI
MEDICAL
PRODUCTS CO., LTD
PROFESSIONAL FACTORY
WITH 20 YEARS' EXPERIENCE

OTHER MEDICAL CONSUMABLES

OTHER MEDICAL CONSUMABLES



03

- 25... COHESIVE BANDAGE
- 26... ADHESIVE BANDAGE
- 27... SPORTS TAPE
- 28... MOSQUITO REPELLENT PATCH
- 28... NASAL INHALANT
- 29... CORN PLASTER
- 30... ESSENTIAL BALM
- 30... ESSENTIAL OIL
- 31... FOAM DRESSING
- 32... SILICONE PLASTER ON POLL
- 33... ONE -PLECE COLOSTOMY BAG
- 34... TWO -PLECE COLOSTOMY BAG



COHESIVE BANDAGE

BL-081

Description

- * Material: Non-woven or cotton
- * Hypo allergenic adhesive, reduce potential skin allergy
- * Adhere to itself, non-sticky to hair or skin
- * Easy to tear, two sides and single-sided adhesive, good elastic.



Item No.	Size	Packing
BL-081	2.5cmx4.5m	24rolls/box, 24boxes/ctn
	5cmx4.5m	12rolls/box, 24boxes/ctn
	7.5cmx4.5m	12rolls/box, 16boxes/ctn
	10cmx4.5m	12rolls/box, 12boxes/ctn
	15cmx4.5m	12rolls/box, 12boxes/ctn



ADHESIVE BANDAGE

BL-081-1

Usage: can bound on any wound part of body

Feature: convenient, popular, non-allergic response





SPORTS TAPE

COTTON SPORTS TAPE (WITHOUT ELASTIC)

(Plain edge/ Zigzag edge)



BL-082



Cotton sport tape is an excellent trainer tape choice. It is extremely soft and conform to different body contours. Strong stickiness, easy hand tear and breathable!


Features:

1. Cotton material
2. Soft and comfortable
3. Good tensile strength
4. Plain edge, easy to tear by hand widthwise
5. Strong and reliable stickiness
6. Constant unwinding tension
7. Allow skin to breathe
8. leaves no residue on body parts

Application:

For fixing muscle when sprains and injuries, also used in fixing ice packs and sports equipment bundle wrapped. It will help control the edema and bleeding.

size	Packing(rolls/boxes/carton)
2.5cmx13.7m	16rolls/box,12boxes/ctn
3.8cmx13.7m	12rolls/box,12boxes/ctn
5.0cmx13.7m	8rolls/box,12boxes/ctn
Custom-made according to customer requirements specification.	



Item No.	Size	Packing
BL-083	2.5cmx5m	24rolls/box,18boxes/ctn
	5cmx5m	9rolls/box,16boxes/ctn
	7.5cmx5m	6rolls/box,24boxes/ctn
	10cmx5m	6rolls/box,18boxes/ctn

SPORTS TAPE(WITH ELASTIC):

BL-83

Kinesiology tape, made of high-quality cotton stretch fabric and medical pressure-sensitive glue, adopts the most advanced "S"shape coating technology in production. It is strong adhesive, difficult falling off, hyposensitivity, soft, air permeability and high ductility. The ductility of the bandage which is considerable equal to the flexibility of muscle is 140%. It is equivalent to an extra muscle on the skin to help muscle contraction, so that is can reduce the burden of muscle without restricting the activities of joints and muscle.

Kinesiology tape is suitable for different joints without drug. Just stick the bandage on the skin when using. By pulling the skin, which administer to circulation of blood and circulation of lymph, balance the uncoordinated muscle. In this way, it can relax the muscle and prevent the muscle pull.



KINESIOLOGY TAPE

BL-84

1. Enhance sports performance
2. Relieve pain
3. Improve circulation
4. Reduce swelling
5. Promote healing
6. Support soft tissue
7. Relax soft tissue
8. Training soft tissue
9. Correct posture
10. Improve action Morphology



MOSQUITO REPELLENT PATCH

BL-085

Paste the patch on the clothes or any location around.
Keep mosquitoes away. Inside or outside your home.
No side effect & effective long lasting up to 12 hours.



Size: $\varnothing 27\text{mm}$, $\varnothing 30\text{mm}$, $\varnothing 38\text{mm}$, $\varnothing 42\text{mm}$ etc.
1Pcs/bag, 2000bags/ctn

NASAL INHALANT

BL-86 N.W. 1g/pcs 1pcs/ blister, 24 blisters/box

Main Ingredients: Menthol, Camphor, etc

Function: Fragrant smell, dispel the wind,
refreshing, freshen brain





CORN PLASTER

BL-087

Reduce friction and pressure on painful calluse, corn, blister and give immediate pain relief.



Item No.	size	Pcs/box	Boxes/middle box	Middle boxes/ctn
BL-086-1	2cm×6cm	5	40	10
BL-086-2	2cm×7cm	4	25	12



ESSENTIAL BALM

BL-088

Main Ingredients: Menthol, Camphor, Peppermint oil, Vaseline etc.

Indications: For refreshment, pain relief, itch relief, Protect against Mosquito bites, Relieve headache caused by heat stroke, and carsick, etc

N.W. : 3g ; 10g ; 19g



Item No.	Size	Packing
BL-088	3g	1pcs/bag, 50 bags/box, 40 boxes/ctn
	3g	1pcs/box, 20 boxes/middle boxes, 40 boxes/ctn
	10g	1pcs/box, 30 boxes/middle boxes, 50 middle boxes/ctn
	19g	1pcs/box, 8 boxes/plastic packing

ESSENTIAL OIL

BL-089

1 bottle/box, 10 boxes/middle boxes

Main Ingredients: Menthol, Camphor, mint oil etc.


Dosage: Put a little on both of the temples or the spot of the skin, 2-3 times per day.

Warning:

For external use only.


Do not apply to wounds or damaged skin.

Avoid contact with eyes.





FOAM DRESSING



BL-097

Number	Size	Packing
01	15*20CM	1pcs/bag,50pcs/box
02	7.5*7.5CM	1pcs/bag,50pcs/box
03	18*18CM	1pcs/bag,50pcs/box
04	15*15CM	1pcs/bag,50pcs/box
05	10*30CM	1pcs/bag,50pcs/box
06	12.5*20CM	1pcs/bag,50pcs/box
07	10*10CM	1pcs/bag,50pcs/box
08	4*5CM	1pcs/bag,50pcs/box

SILICONE PLASTER ON ROLL




BL-098

Number	size	Packing
01	7.5*4CM	one roll one box,50boxes/ctn
02	4*15CM	one roll one box,50boxes/ctn
03	3.5CM*1.5M	one roll one box,50boxes/ctn
04	4CM*3M	one roll one box,50boxes/ctn



THE LAYERS CAN BE CUSTOMIZED



ONE-PIECE PRODUCT ANALYSIS PROFILE

1. Release film (Strong adhesive, easy to peel off)
2. Medical hydrocolloid (Using medical hydrocolloid materials, good stickiness, long-lasting and skin-friendly, no residue on the skin)
3. Spunlace extender (Reinforce the baseplate, prevent warping, not prone to allergies)
4. Medical non-woven fabric (Soft, comfortable and ventilate, reduce friction noise, avoid embarrassment)
5. EVOH film (7-layer film co-extruded, strong and durable)
6. Back window (Double layer medical non-woven fabric, reduce friction sound, more concealed)
7. Exhaust valve (Quickly discharge the flatulence in the bag)
8. Carbon filter (Good adsorption, effectively filter excrement odor)
9. Sealing method (The sealing operation is convenient and the sealing is tight)

SEALING METHOD

There are 5 types of sealing methods for ostomy bags, 4 types for colostomy bag, namely: clamp, hook and loop, Velcro and medical twist tie; 1 type for urostomy bag: drain valve.



Clamp Hook and loop Velcro Medical twist tie Drain valve

COLOSTOMY CARE

Ostomy bag is used to collect stool and urine from an ostomy patient. There are two types, single-piece bags and two-piece bags according to their structure. A one-piece ostomy bag is used for various ostomies (colon, ileum, urinary tract). It is easy to use. The two-piece bag has two parts, Skin Barrier and Colostomy pouch. It is to avoid damage to the stoma and surrounding skin, and only need to replace the bag body.

Colostomy pouch
Size: 480ml, 310ml, 210ml, 450ml

BL96



Skin Barrier



Ostomy Attachment





ONE- PIECE COLOSTOMY BAG

9501



9502



9503



9504



9505



9506



BL-95

TWO- PIECE COLOSTOMY BAG

9601



9602



MODEL	FLANGE	CAN BE MANUFACTURED	SIMPLE PACK
9601	45	45 57 70 FLANGE	20PCS/BOX,40BOXES/CARTON
9602	57	45 57 70 FLANGE	20PCS/BOX,40BOXES/CARTON

9603



9604



MODEL	FLANGE	CAN BE MANUFACTURED	SIMPLE PACK
9603	70	45 57 70 FLANGE	20PCS/BOX,40BOXES/CARTON
9604	80	45 57 70 FLANGE	20PCS/BOX,40BOXES/CARTON

BL-96